




-  Bedside blood gas analysis
-  Lightweight and portable
-  22 parameters in one test



Blood Gas Analyzer

EG-i Series

EAGLENOS Co., Ltd.

Blood Gas Analyzer



Accurate, reliable and maintenance-free

Self-calibration for each test
Accurate results in approximately 5 minutes
Dry chemistry method, no reagent pack required, no carry-over contamination



Lightweight and portable

Rechargeable lithium battery: more than 50 tests
Size: 235mm×210mm×160mm
Weight: 3±0.5 kg (including battery)



Easy-to-use

Quick start tutorial
8-inch full HD touchscreen



Smart identification in cartridge

Feedback of cartridge insertion
Identification of the cartridge expiration date



Automatic quality control

Regular QC reminders
Power on Self-Test
Dual quality control: electronic simulator and controls



Blood Gas Test Cartridge 25 tests/kit



22 parameters of blood gases, electrolytes and metabolites

10 measured parameters: K^+ , Na^+ , Cl^- , iCa^{2+} , pH, pCO_2 , pO_2 , Glu, Lac, Hct
12 calculated parameters: HCO_3^- , TCO_2 , BE-ecf, SO_2 , Gap, Gap- K^+ , $pO_2(T)$, $pCO_2(T)$, pH(T), BE(b), $pO_2(A-a)$, $pO_2(a/A)$



Precise industrial design

High performance sensors
Microfluidic technology and integrated flow control system
Unique sample quantification design to precisely control sample volume



Instant detection, cost-effective and environment friendly

Rapid analysis of blood (venous/arterial) without sample processing
Single-use cartridge to avoid waste and cross-contamination



Accurate

Automatic calibration upon each test to ensure accurate results



Test parameters and clinical significance

Electrolytes

Potassium ion (K^+)

Even small changes in extracellular K^+ concentration will have significant effects on the transmembrane potential gradient, and thereby the function of neuromuscular and cardiac tissues.

Sodium ion (Na^+)

As the most abundant extracellular fluid solute, Na^+ is the major determinant of its osmolality and thereby the principal determinant of water distribution between the intracellular and extracellular compartments. This highlights the role of Na^+ in the maintenance of blood volume and thereby blood pressure.

Chloride ion (Cl^-)

As the second most abundant extracellular fluid ion after Na^+ , and the most abundant extracellular fluid anion, Cl^- is essential for the maintenance of normal plasma osmolality.

Free calcium ion (iCa^{2+})

The maintenance of iCa^{2+} within normal limits is not only important for the structural integrity of bones but for a range of physiological functions, including: hemostasis, cardiac and skeletal muscle cell contraction, neuromuscular transmission and action of many hormones (calcium-signaling).

pH, Blood Gas

Acidity and alkalinity (pH)

The pH level is an indicator of the acidity and alkalinity of the blood. The abnormal pH level indicates an acid-base imbalance.

Carbon dioxide partial pressure (pCO_2)

pCO_2 is the partial pressure generated by physically dissolved CO_2 molecules in the blood and is an important indicator of the effectiveness of alveolar ventilation.

Oxygen partial pressure (pO_2)

pO_2 is the partial pressure generated by the physically dissolved O_2 molecules in the blood and reflects the oxygen uptake by the pulmonary capillary blood.

Biochemical metabolites/Hematocrit

Glucose concentration (Glu)

Glucose is the primary source of energy for the organism and the only source of specialized nutrition for brain tissue. The measurement of blood glucose levels is extremely important for the diagnosis and treatment of patients with diabetes and hypoglycemia.

Lactic acid concentration (Lac)

Lactate is an indicator to assess the degree of tissue hypoperfusion and cellular hypoxia.

Hematocrit (Hct)

The percentage of red blood cells to whole blood volume is the main indicator of blood viscosity, anemia, severe blood loss and the body's ability to transfuse oxygen.

Application Scenario



Emergency Department

Monitoring of blood gases, metabolites and electrolyte parameters for emergency critical care patients such as poisoning, coma and convulsions



ICU

Monitoring of blood gases, metabolites, electrolytes and acid-base balance for critically ill patients



Anesthesiology

Blood gases and electrolyte monitoring during surgical anesthesia (preoperative, intraoperative, postoperative)



Surgery

Blood gases, electrolyte and acid-base balance monitoring during surgery



Pediatrics/ Fever Clinic

Electrolyte monitoring for outpatient infusion



Respiratory

Determine the classification of respiratory failure, assessment of the degree of hypoxia, and guidance on ventilator adjustment



Gastroenterology

Electrolyte and acid-base balance monitoring in patients with diarrhea and vomiting with coma



Cardiology

Blood gases and electrolyte monitoring in patients with heart attack



Endocrinology

Blood gases, metabolites, electrolyte and acid-base balance monitoring in patients with diabetic ketoacidosis



Neurology

Blood gases and acid-base balance monitoring in patients with stroke and coma



Obstetrics and Gynecology/ Neonatology

Diagnosis of neonatal asphyxia, assessment of neonatal pulmonary disease and monitoring of acid-base balance


Specification

Parameter	Accuracy	Range	Precision
K⁺	±3.0%	2.0-9.0 mmol/L	≤1.5%
Na⁺	±3.0%	100-180 mmol/L	≤1.5%
Cl⁻	±3.0%	65-140 mmol/L	≤1.5%
iCa²⁺	≥1.00 mmol/L, ±5.0%; <1.00 mmol/L, ±0.05 mmol/L	0.25-2.50 mmol/L	≤1.5%
pH	±0.04	6.500-8.000	≤0.02
pCO₂	≥62.0mmHg, ±8.0%; <62.0mmHg, ±5.0 mmHg	10.0-150.0 mmHg	≥62.0mmHg, CV≤4.0%; <62.0mmHg, SD≤2.5 mmHg
pO₂	≥50mmHg, ±15.0%; <50mmHg, ±7.5 mmHg	10-700 mmHg	≤5.0%
Glu	≥4.0 mmol/L, ±10.0%; <4.0 mmol/L, ±0.4 mmol/L	1.1-38.0 mmol/L	≥4.0 mmol/L, CV≤5.0%; <4.0 mmol/L, SD≤0.2 mmol/L
Lac	≥5.00 mmol/L, ±12.0%; <5.00 mmol/L, ±0.6 mmol/L	0.50-20.00 mmol/L	≥5.00 mmol/L, CV≤6.0%; <5.00 mmol/L, SD≤0.3 mmol/L
Hct	≥50%PCV, ±6.0%; <50%PCV, ±3%PCV	10%-70%PCV	≥50%PCV, CV≤3%; <50%PCV, SD≤1.5%PCV

Operation instruction



Step1: Scan the barcode on the cartridge pouch.



Step2: Fill the inlet with sample and slide the cap to seal the sample inlet.



Step3: Insert the cartridge until it clicks. Wait for the results.

Company info

Nanjing Eaglenos Co., Ltd.

Tel: +86-025-400-019-0069

Address: B2-2, Treehouse Headquarters Cluster, No.73 Tanmi Road Jiangbei New Area, Nanjing, Jiangsu 210044, China

Website: www.eaglenos.com

If you need to know the intended use of the product, precautions and contraindications, please refer to the instructions.



*For health care professionals only