


CBS-30 Electrolyte Analyzer (K^+ Na^+ Cl^-)



Products Characteristics

- Advanced Multi- Biosensor Technology (MBT) permits maintenance-free electrode.
- Close reagent pack is easy to replace.
- Multi-level calibrators and controls secure results accuracy.
- User friendly operation and quick trouble shooting.
- High throughput gives each test result within half a minute.
- Up to 6000 sample results storage and easy to review.
- 5.6 inch touch screen makes operation more convenient

Specifications

Method	Multi-Biosensor Technology			
Sample Size	60-120 μ l			
Speed	≥ 120 tests/hour			
Measure scope & Precision	Performance	Detection Range(mmol/L)	Precision(CV)	Resolution(mmol/L)
	K ⁺	0.50 ~ 20.00	$\leq 1.0\%$	0.01
	Na ⁺	15.0 ~ 200.0	$\leq 1.0\%$	0.1
	Cl ⁻	15.0 ~ 200.0	$\leq 1.0\%$	0.1
Printer	Inner thermal printer			
Expand Function	RS232 & USB Connector, Computer and Barcode reader available			
Ambient Conditions	0-40 $^{\circ}$ C , Relative humidity $\leq 85\%$			
Power	AC 80V ~ 260V , 50-60Hz, ≤ 80 W			
Display	Color display with 5.6 inch touch screen			



B&E Diagnostics Inc.

700 Nicholas Blvd, Ste 204 Elk Grove Village, IL 60007, USA

Tel: +1 630 974 5159 Fax: +1 630 974 5160

info@besic.com www.besic.com