

Specification

Dimensions(mm)	260, 240, 140
Weight	2.6 Kg
Data Storage	8000 tests result
Power Adapter	AC 100~240V, 50/60 Hz
Data Output	Onboard screen/Printer/PC/LIS
Rated Power	36W
Displayer	8 inch high resolution touch screen
QC Code Reading	RFID
Repeatability	CV<3%
Accuracy	± 3%
Stability	CV<3%
Excitataion Light Source	365 nm
Detector	615 nm
Humidity	≤80%




Lamuno Plus

Multi Scenes Usage


Fluorescence Immunoassay Analyzer



Highlight Feature



Testing Process



Application scenarios



Tests List

Cardiac Marker	Tumor Marker	Hormone	Inflammation Detection
NT-proBNP	Ferritin	β -HCG	Hs-CRP + CRP
CK-MB/CTnI/Myo	PSA	LH	PCT
cTnI	AFP	FSH	SAA
cTnT	f-PSA	PRL	IL-6
Myo	CEA	E2	
CK-MB	CA125	AMH	
D-Dimer	CA153	Tes	
sST2	CA199	Prog	
			COVID-19
			COVID19 Neutralizing Antibody
			COVID19 Antigen
Digestive Function	Thyroid Function	Renal Function	Others
G17	T3	NGAL	IgE
Cal	T4	MAU	S100- β
FOB	FT3		HbA1c
PGI/PGII	FT4		25-OH-VD
	Cor		
	TSH		
		Infectious	
		HCV	
		HBsAg	